

Custom Recombinant Protein Services

Comprehensive Clone to Protein Service

OriGene offers comprehensive, high-quality recombinant custom protein services from different expression systems, in addition to the standard catalog items.

- ✓ Bulk quantity (µg to mg)
- ✓ Various hosts - HEK293, CHO, and *E. coli*
- ✓ Custom Tags or Tag-less
- ✓ Mutation, buffers
- ✓ SEC testing
- ✓ Custom concentrations
- ✓ Dedicated PhD level project manager




All our proteins are proudly produced in the USA

2,000+
Successful
Projects


>500
Satisfied
Customers

25+
Years
Experience


Successful Examples



Biotinylated Protein



Mutant Protein



Tag-Free Protein



Fluorescent-labeled Protein



OriGene Headquarters
9620 Medical Center Dr
Rockville, MD 20850
+1-301-340-3188
custsupport@origene.com
techsupport@origene.com
sales@origene.com

OriGene EU
SchillerstraBe 5
32052 Herford, Germany
+49-5221-346060-0
info-de@origene.com

www.origene.com

ORIGENE
EMPOWER YOUR RESEARCH