


cDNA clone vectors

TrueClone vectors (non-tagged)

Most of the non-tagged cDNA clones were cloned in pCMV6-XL4, 5 and 6 vectors. The cDNA inserts are easily shuttled into pCMV6-Neo vector which contains a neomycin mammalian selection marker.



TrueORF Vectors (tagged)

	pCMV6-Entry	pCMV6-AC-GFP	pLenti-C-Myc-DDK	pLenti-C-mGFP
Vector SKU	PS100001	PS100010	PS100064	PS100071
Vector Map				
Fusion Tag	C-terminal Myc-DDK	C-terminal tGFP	C-terminal Myc-DDK	C-terminal mGFP
Ideal For	Detection and purification of over-expressed protein	Tracking the over-expressed protein in transfected cells	In hard-to-transfect cells: Detection and purification of over-expressed protein	Tracking the over-expressed protein in transfected cells
ORF Clone SKUs	RC2XXXXX	RG2XXXXX	RC2XXXXX1	RC2XXXXX2