


OriGene mAb Services

Key Services

- Rabbit mAb via B-cell cloning
- Mouse hybridoma development
- Specificity test using 10K chip




Features:

- ✓ Proven track record – Successfully produced >20,000 mouse mAbs & >100 rabbit mAbs
- ✓ Guaranteed lead time - Mouse mAb (6~10 months) & rabbit mAb (4~6 months)
- ✓ Proprietary platform for Ab specificity test High-density protein chip (10,000)

Mouse mAb


High sensitivity for native epitopes



Positive staining of human liver carcinoma cells by PD-L1 mouse mAb (Cat#UM800120)

Rabbit mAb

A great tool for mouse model




Positive staining of mouse APC immune cells by PD-L1 rabbit mAb (Cat#TA591004)

For more information, please visit

<https://www.origene.com/services/antibody-development>

Recombinant Rabbit mAb

❖ Rabbit B cell cloning and recombinant antibody production



Advantages of recombinant rabbit mAbs:

- Accelerated development time
- Higher sensitivity for mutants or phosphorylated proteins
- Enhanced reproducibility than hybridoma

For inquiry and order, please contact sales@origene.com