



B-cell Lineage Markers




Popular B-cell Markers

CD19




IHC staining of carcinoma of human lung tissue using anti-CD19 mouse monoclonal antibody

CD10



FC analysis of living 293T cells transfected with MME overexpression plasmid using anti-MME


SDC1



IHC staining of human appendix tissue within the normal limits using anti-SDC1

TA506236S


CD20



IHC staining of human tonsil using anti-CD20

UM870128


CD21



IHC staining of human lymphoma using anti-CD21

TA813767S

CD22



IHC staining of human lymphoma using anti-CD22

TA800385S

TA814111S

TA506412S

B-cell Lineage Marker Antibody Panel

Biomarker	SDC1	MS4A1	MME	CR2	CD22	CD19
Catalog-No.	TA813767S	TA800385S	UM870128	TA814111S	TA506412S	TA506236S
Description	SDC1 mouse monoclonal antibody, clone OTI3G1	MS4A1 (CD20) mouse monoclonal antibody, clone OTI4B4	MME/CD10 mouse monoclonal antibody, clone UMAB236	CD21 mouse monoclonal antibody, clone OTI9D9	CD22 mouse monoclonal antibody, clone OTI4C3	CD19 mouse monoclonal antibody, clone OTI3B10

Get a sampler antibody Kit for B-cell lineage marker, contact our tech support team at techsupport@origene.com

B-Cell Markers

Target	Reactivity	Application	SKU
CD5	Human	FC,IHC, WB	TA802427S
CD10	Human, Mouse, Rat	10k-ChIP, FC, IHC, WB	UM870128
CD19	Human	FC,IF, IHC, WB	TA506236S
CD21	Human	IHC	TA814111S
CD22	Human, Dog, Rat, Monkey, Mouse	IHC, WB	TA506412S
CD23	Human	IF, IHC, WB	UM870031
CD24	Human	IHC	TA366996
CD27	Human, Rat, Mouse	IHC, WB	TA806591S
CD34	Human	IHC, WB	TA808885S
CD38	Human	IHC, WB	TA804538S
CD40	Human	IHC, WB	TA807569S
CD53	Human	FC, IF, IHC, IP, WB	TA386337
CD69	Human	FC, FN, IHC, WB	AM03132PU-N
CD72	Human	FC	SM115PE
CD80	Human	IF, IHC, WB	TA501576S
CD86	Human	FC, WB	TA813079S
CD93	Human	IHC	TA322420
CD95	Human	WB	TA506937S



For more details on
markers visit us at
www.origene.com