

siTran 1.0 Citations

Update – 02/2013

Taxifolin Enhances Andrographolide-Induced Mitotic Arrest and Apoptosis in Human Prostate Cancer Cells via Spindle Assembly Checkpoint Activation

Zhang ZR, Al Zaharna M, Wong MM, Chiu SK, Cheung HY. PLoS One. 2013;8(1):e54577

[Successful transfection of DU145 cells]

Inhibition of the de-myelinating properties of Aicardi-Goutières Syndrome lymphocytes by cathepsin D silencing

Pulliero A, Marengo B, Longobardi M, Fazzi E, Orcesi S, Olivieri I, Cereda C, Domenicotti C, Balottin U, Izzotti A. Biochem Biophys Res Commun. 2013 Jan 18;430(3):957-62

[Successful transfection of lymphocytes]

A possible gene silencing mechanism: Hypermethylation of the Keap1 promoter abrogates binding of the transcription factor Sp1 in lung cancer cells

D Guo, B Wu, J Yan, X Li, H Sun, D Zhou Biochem Biophys Res Commun. 2012 Nov 9;428(1):80-5

[Successful transfection of BEAS 2B cells]

Human β -defensin-2 enhances IFN- and IL-10 production and suppresses IL-17 production in T cells

Naoko Kanda, Masahiro Kamata, Yayoi Tada, Takeko Ishikawa, Shinichi Sato, and Shinichi Watanabe J. Leukoc. Biol., Jun 2011; 89: 935 – 944

[Successful transfection of T-cells]