

**PRODUCT: Nuclease Free Duplex Buffer (2 ml)**

Catalogue Number: SR30005

This buffer is recommended for reconstitution of OriGene's lyophilized Trilencer siRNA products and the controls.

## Product specification:

|                               | <b>Specification</b> |
|-------------------------------|----------------------|
| Visual Description            | Colorless Solution   |
| Determination of Particulates | No Particulates      |
| Determination of pH           | 7.3-7.8              |
| Endotoxin Contamination       | <0.25EU/mL           |
| RNaseA                        | <1 ng/mL             |
| RNaseT1                       | <12.5 ng/mL          |
| RNase1                        | <70 ng/mL            |

Storage condition: -20C.