

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Anti-Fatty Acid Binding Protein (FABP)
-HRP Conjugate

Catalog Number: 80099

Lot Number:

Quantity: _____ ml

Antibody Concentration ($A_{280\text{nm}}$): _____ mg/ml

Product Form: In 0.015 M KPO_4 Buffer, 0.85 NaCl, pH = 7.20,
containing proteins and preservatives.

Storage: At 2-8° C. DO NOT FREEZE. The conjugate may
become inactivated upon freezing.

Pairing Information: This conjugate pairs with goat polyclonal anti-fatty acid
binding protein (P/N: 70385) in the human Fatty Acid
Binding Protein sandwich immunoassay.