

Product Code: B9085MA06-MA
Product Name: Anti-Human Cardiac Troponin I (24-40)
Description: Troponin I Monoclonal, Protein A purification
Lot Number:
Site of Manufacture: Newark, DE

Antibody Purification

Biological Source: Mouse, ascites was obtained from healthy animals under NIH guidelines (*In vivo*)
Product Form: Liquid
Buffer: 20mM Sodium phosphate, 150mM Sodium chloride, pH 7.2
Specificity: Human Cardiac Troponin I (24-40)
Clone: TPC 110
Country of Origin: USA
Filtration: 0.2µm
Date of Filtration:
Purification: Protein A
Storage Conditions: 4°C (Range 2°C - 8°C)

Test:	Method:	Document:	Specification:	Results:
Protein Concentration	A280	G3476	3-5 mg/mL	
Purity	HPLC	G3331	≥ 95%	≥ 95%
Isoelectric Point	Isoelectric Focusing	G3305	Match Reference	6.3 – 6.7 pI range (matches reference 700.389-3-5)
Isotype	ISOSTrip	G3346	Mouse IgG1, Kappa	Mouse IgG1, kappa

Intended Use: For Further Manufacturing or Research Use Only. Product meets European Union requirements for treated technical blood products. Not for Human or Animal Consumption.

Prepared By: _____

Reviewed By: _____

Date: _____

Date: _____