

Product Code: D4160MA01-MA
Description: Troponin I Monoclonal (41-60), Protein A Purification
Lot Number:
Site of Manufacture: Newark, DE

Antibody Purification

Biological Source: Mouse, USA. Hybridoma Cell Culture Supernatant
Product Form: Liquid in 20mM Sodium phosphate, 150mM Sodium chloride, pH 7.2
Specificity: Human Cardiac Troponin I (41-60)
Clone: 10C9
Filtration: 0.2µm
Date of Filtration:
Storage Conditions: 4°C (Range 2°C - 8°C)

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	HPLC	G3331	≥ 95%
Isoelectric Point	Isoelectric Focusing	G3305	5.9 – 6.4 pI range
Isotype	ISOStrip™	G3346	Mouse IgG1, kappa

Prepared By: _____

Date Reviewed: _____

Reviewed By: _____

Date Reviewed: _____