

Product Name: NT-proBNP mouse monoclonal antibody ,clone 2F4
Product: Purified Monoclonal Antibody **Catalog No.:** D601
Lot Number: **Concentration:**
EXP Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

| Test: | Method: | Document: | Results: |
|-----------------------|--------------------|-----------|------------|
| Protein Concentration | BCA | G3576 | |
| Purity | SDS-Page (reduced) | G3302 | > 95% |
| Isotype | ISOSTrip™ | G3346 | Mouse IgG1 |
| Application | WB | Datasheet | PASS |

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date