



Certificate of Analysis

Product Name: NT-proBNP mouse monoclonal antibody, clone 10A2
Product: Purified Monoclonal Antibody **Catalog No.:** D606
Lot Number: **Concentration:**
EXP Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: Freeze-dried in PBS, BSA/glycerol-free, reconstitute with PBS (pH7.3)
Storage Conditions: -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3476	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	ELISA	Datasheet	

Prepared By:

Reviewed By:

9620 Medical Center Dr. | Rockville, MD20850 USA | +1 (888) 267-4436 | +1 (301) 340-3188

www.origene.com custsupport@origene.com