

**Product Name:** cTnI-biotin  
**Product:** Biotin-coupled Monoclonal Antibody  
**Lot Number:**  
**Production Date:**  
**Exp Date:**

**Catalog No.:** DM0022-B14  
**Concentration:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Application	ELISA	Datasheet	

**Prepared By:** \_\_\_\_\_  
**Quality Control-Date**

**Reviewed By:** \_\_\_\_\_  
**Quality Assurance-Date**