

Product Name: Troponin I (24-40) Rabbit monoclonal antibody, clone OTIR1D6
Product: Purified Monoclonal Antibody **Catalog No.:** DM0022-S1
Lot Number: **Concentration:**
EXP Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein A Sepharose
Buffer/Components: PBS (pH7.0), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOStrip™	G3346	
Application	Elisa	Datasheet	

Prepared By:**Reviewed By:**

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.comcustsupport@origene.com