

Product Name: BNP Rabbit monoclonal antibody, clone OTIR2A6
Product: Purified Monoclonal Antibody **Catalog No.:** DM0025-S2
Lot Number: **Concentration:**
EXP Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein A Sepharose
Buffer/Components: PBS (pH7.0), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

| Test: | Method: | Document: | Results: |
|-----------------------|--------------------|-----------|----------|
| Protein Concentration | BCA | G3576 | |
| Purity | SDS-Page (reduced) | G3302 | |
| Isotype | ISOSTrip™ | G3346 | |
| Application | Elisa | Datasheet | |

Prepared By:**Reviewed By:**

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.comcustsupport@origene.com