

Product Name: ST2 mouse monoclonal antibody ,clone 21A7
Product: Purified Monoclonal Antibody
Lot Number: **Catalog No.:** TA810124
EXP Date: **Concentration:**

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	> 95%
Isotype	ISOSTrip™	G3346	Mouse IgG1
Application	WB	Datasheet	PASS

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date