

Product Name: SELENBP1
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: TP302108
Concentration :
EXP date:

Antibody Purification

Molecular Weight: 52.2 kDa
Expression Host: HEK293T
Starting Material: Cell Lysate
Buffer/Components: 25mM Tris-HCl, pH7.3, 100mM glycine, 10% glycerol.
Storage Conditions: -80°C, avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA methods		
Purity	SDS-Page(reduced)	G3302	

Certified by: _____ **Date:** _____