

Certificate of Analysis

Product Information

Product name:NT-proBNP protein

Catalog No:TP750122

Lot Number:

Quantity:

Concentration:

Buffer Formulation:PBS,pH7.4,10% glycerol

Storage Conditions: Store at -80°C. Avoid repeated freeze-thaw cycles. Stable for at least 3 years from receipt of products under proper storage and handling conditions.

Quality Control

Test Items	Specifications	Results:
Protein Concentration	Microplate BCA Protein Assay	Conformed
Purity(SDS-PAGE)	> ____%	Conformed

Certified by: _____

Date: _____