

## Certificate of Analysis

---

### Product Information

Product name:CKMB protein

Catalog No:TP750145

Lot Number:

Quantity:

Concentration:

Buffer Formulation:10 mM Tris, pH 7.2, 5 mM mercaptoethanol, 50% glycerol.

Storage Conditions: Store at -80°C. Avoid repeated freeze-thaw cycles. Stable for at least 3 years from receipt of products under proper storage and handling conditions.

---

### Quality Control

Test Items	Specifications	Results:
Protein Concentration	Microplate BCA Protein Assay	Conformed
Purity(SDS-PAGE)	> ____%	Conformed

Certified by: \_\_\_\_\_

Date: \_\_\_\_\_