

Product Name: PCNA
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: TP750150
Concentration :
EXP date:

Antibody Purification

Molecular Weight: 28.8KDa
Expression Host: E. coli
Starting Material: Lysate
Buffer/Components: 20mM Tris, pH7.5, 2mM EDTA, 20% Glycerol.
Storage Conditions: -80°C, avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA methods		
Purity	SDS-Page(reduced)	G3302	

Certified by: _____ **Date:** _____

No.168 Mashang Meiliang Road | Wuxi , Jiangsu province, China | 0510-81830666