

## Certificate of Analysis

---

### Product Information

Product name:cTnT protein

Catalog No:TP750172

Lot Number:

Quantity:

Concentration:

Buffer Formulation:Lyophilized from 0.01 M HCl,reconstitute with Tris/urea buffer (20 mM Tris pH7.5, 7 M urea, 5 mM EDTA, 15 mM  $\beta$ -mercaptoethanol).

Storage Conditions: Store at -80°C. Avoid repeated freeze-thaw cycles. Stable for at least 3 years from receipt of products under proper storage and handling conditions.

---

### Quality Control

Test Items	Specifications	Results:
Protein Concentration	Microplate BCA Protein Assay	Conformed
Purity(SDS-PAGE)	> ____%	Conformed

Certified by:\_\_\_\_\_

Date:\_\_\_\_\_

---