



Certificate of Analysis

Product Name: Troponin (I-C) Protein
Catalog No.: TP810062
Lot Number:
EXP Date:

Concentration:

Antibody Purification

Buffer/Components: Tris Buffer, pH 8.0, containing 0.45% NaCl and 50% glycerol
Storage Conditions: -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Application	ELISA	Datasheet	

Prepared By:

Reviewed By:

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.com

custsupport@origene.com