

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Polyclonal Anti-Human Prostate Specific Antigen Total PSA (Total IgG Fractions)

Catalog Number: 70339 (1110)

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280} nm): _____ mg/ml

Product Form: DEAE column chromatography purified immunoglobulin fraction in 0.015M KPO_4 buffer, pH = 7.40, containing 0.85% NaCl, 0.1% NaN_3 .

Specificity: Human Prostate-specific antigen

Cross-Reactivity: No cross-reactivity with Human Albumin, AFP, CEA, Ferritin, hCG, PAP, CA 125, CA 15-3, or CA 19-9.

Storage: At 2-8° C.

Suggested Use: For solid phase coating (capture).

Pairing Information: This antibody pairs with mouse monoclonal anti-human PSA (P/N: 70338) in the human Total PSA sandwich immunoassay.