

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	<b>Rabbit Polyclonal Anti Human Ferritin, Affinity Purified</b>
<b>Catalog Number:</b>	70643
<b>Lot Number:</b>	
<b>Quantity:</b>	_____ mg (= _____ ml)
<b>Antibody Concentration (<math>A_{280nm}</math>):</b>	_____ mg/ml
<b>Product Form:</b>	Ferritin affinity column chromatography purified antibody in 0.015 M potassium phosphate buffer, pH=7.40 containing 0.85% NaCl, and 0.1% NaN <sub>3</sub> .
<b>Storage:</b>	At 2-8°C.
<b>Specificity:</b>	This antibody reacts with human ferritin (liver and spleen).
<b>Suggested Use:</b>	For enzyme conjugation (detection) or solid phase coating (capture).
<b>Pairing Information:</b>	This antibody pairs with the following mouse monoclonal antibodies in the human Ferritin Sandwich Immunoassay. <ol style="list-style-type: none"><li>1. Mouse monoclonal anti-ferritin, P/N: 70226 (capture or detection).</li><li>2. Mouse monoclonal anti-ferritin, P/N: 70641 (capture).</li></ol>