

**Product Name:** PSMA mouse monoclonal antibody ,clone 2E1  
**Product:** Purified Monoclonal Antibody **Catalog No.:** D121  
**Lot Number:** **Concentration:**  
**EXP Date:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	> 95%
Isotype	ISOSTrip™	G3346	Mouse IgG2b
Application	ELISA	Datasheet	PASS

Prepared By: \_\_\_\_\_  
Quality Control-Date

Reviewed By: \_\_\_\_\_  
Quality Assurance-Date