

**Product Name:** cyfra21-l mouse monoclonal antibody ,clone 9A9  
**Product:** Purified Monoclonal Antibody **Catalog No.:** DM0019-B1  
**EXP Date:** May.2023 **Concentration:** 3.48mg/ml

### **Antibody Purification**

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C, Do Not freeze.

<b>Test:</b>	<b>Method:</b>	<b>Document:</b>	<b>Results:</b>
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	WB	Datasheet	
Application	IHC	Datasheet	

**Prepared By:** \_\_\_\_\_  
Quality Control-Date

**Reviewed By:** \_\_\_\_\_  
Quality Assurance-Date