

Product Name: cyfra21-l mouse monoclonal antibody ,clone 5H6
Product: Purified Monoclonal Antibody **Catalog No.:** DM0019-B2
EXP Date: June.2023 **Concentration:** 5.89mg/ml
Lot Number:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	WB	Datasheet	
Application	IHC	Datasheet	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date