



Certificate of Analysis

Product Name: Mouse Anti-human Cyfra21-1 monoclonal antibody ,clone 5C10
Product: Purified Monoclonal Antibody **Catalog No.:** DM0019-B5
Lot Number: V0102 **Concentration:** 4.58mg/ml
EXP Date: July.2027

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	4.58mg/ml
Purity	SDS-Page (reduced)	G3302	> 95%
Isotype	ISOSTrip™	G3346	Mouse IgG1
Application	ELISA	Datasheet	PASS

Prepared By: Tim Li 2022.07.15
Quality Control-Date

Reviewed By: Lixiaom 2022.07.22
Quality Assurance-Date

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.com

custsupport@origene.com