

Product Name: Mouse Anti-human GPC3 monoclonal antibody, clone 2H10
Product: Purified Monoclonal Antibody **Catalog No.:** DM0064-B2
Lot Number: **Concentration:**
Production Date:
EXP Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN₃
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	Elisa	Datasheet	

Prepared By:**Reviewed By:**

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.comcustsupport@origene.com