

Product Name: G-17 Protein
Product: Purified Monoclonal Antibody
Lot Number:
EXP Date:

Catalog No.: DP0004-P1
Concentration:

Antibody Purification

Buffer/Components: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Conditions: -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	> 95%
Application	ELISA	Datasheet	PASS

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date