

**Product Name:** ProGRP  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Manufactured:**

**Catalog No.:** DP0016-P1  
**Concentration :**  
**EXP date:**

## Antibody Purification

**Buffer/Components:** 20mM PB, 75mM NaCl, 0.05% NaN<sub>3</sub>, 50% glycerol  
**Storage Conditions:** -20 °C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA methods		
Purity	SDS-Page(reduced)	G3302	

**Certified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

No.168 Mashang Meiliang Road | Wuxi , Jiangsu province, China | 0510-81830666

[www.origene.com](http://www.origene.com)

[custsupport@origene.com](mailto:custsupport@origene.com)