



Certificate of Analysis

Product Name: CA242 protein
Catalog No.: DP0022-P1
Lot Number:
EXP Date:

Concentration:

Antibody Purification

Buffer/Components: PBS (pH7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage Conditions: -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Application	ELISA	Datasheet	

*The data from the raw material without adding the dilution buffer.

Prepared By:

Reviewed By:

No.168 Mashan Mellang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.com

custsupport@origene.com