



## Certificate of Analysis

**Product Name:** CA153 protein  
**Catalog No.:** DP0034-P1  
**Lot Number:**  
**EXP Date:**

**Concentration:**

### Antibody Purification

**Buffer/Components:** PBS (pH7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.  
**Storage Conditions:** -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Application	ELISA	Datasheet	

**Prepared By:**

**Reviewed By:**

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

[www.origene.com](http://www.origene.com)

[custsupport@origene.com](mailto:custsupport@origene.com)