

TICertificate of Analysis

Product Information

Product name: GOLM1 protein

Catalog No: TP314745

Lot Number:

Quantity:

Concentration:

Buffer Formulation: 25 mM Tris-HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Storage Conditions: Store at - 80 °C. Avoid repeated freeze-thaw cycles. Stable for at least 3 years from receipt of products under proper storage and handling conditions.

Quality Control

Test Items	Specifications	Results:
Protein Concentration	Microplate BCA Protein Assay	Conformed
Purity(SDS-PAGE)	> ___%	Conformed

Certified by: _____

Date: _____