

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Human CRP, Affinity Purified.

Code: 70189

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: CRP antigen affinity column chromatography purified antibodies in 15 mM potassium phosphate buffer, pH = 7.40, containing 0.85% NaCl and 0.1% NaN₃.

Storage: 2-8° C

Specificity: This antibody reacts with human CRP.

Suggested Use: For enzyme conjugation (detection).

Pairing Information: This antibody pairs with goat anti-CRP (P/N: 70243) in the human CRP sandwich immunoassay.