

**Product Name:** IL17A mouse monoclonal antibody, clone 1G8  
**Product:** Purified Monoclonal Antibody      **Catalog No.:** DM0032-B2  
**Lot Number:**      **Concentration:**  
**EXP Date:**

**Antibody Purification**

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

<b>Test:</b>	<b>Method:</b>	<b>Document:</b>	<b>Results:</b>
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	Elisa	Datasheet	

**Prepared By:****Reviewed By:**

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666