

## TICertificate of Analysis

---

### Product Information

Product name: CRP

Catalog No: DP0002-N1

Lot Number: BA10013

Quantity: 50 ug

Concentration: 2.85ug/ul

Buffer Formulation: 25 mM Tris.HCl, pH 7.3, 100 mM glycine

Tags: Nature protein purified from human serum

Storage Conditions: Store at -20°C. Avoid repeated freeze-thaw cycles. Stable for at least 3 months from receipt of products under proper storage and handling conditions.

---

### Quality Control

Test Items	Specifications	Results:
Protein Concentration	Microplate BCA Protein Assay	Conformed
Purity(SDS-PAGE)	>80%	Conformed

Certified by: \_\_\_\_\_

Date: \_\_\_\_\_