

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	<b>Goat Polyclonal Anti-Human Whole hCG, Affinity Purified</b>
<b>Catalog Number:</b>	70396
<b>Lot Number:</b>	
<b>Quantity:</b>	_____ mg (= _____ ml)
<b>Antibody Concentration (<math>A_{280nm}</math>):</b>	_____ mg/ml
<b>Product Form:</b>	In 0.015 M potassium phosphate buffer, 0.85% NaCl, 0.1% $NaN_3$ , pH = 7.40
<b>Purity:</b>	Greater than 98% pure by PAGE
<b>Storage:</b>	At 2-8°C.
<b>Specificity:</b>	This antibody reacts with alpha chain and beta chain of the human whole hCG molecule.
<b>Pairing Information:</b>	<ol style="list-style-type: none"><li>1. This antibody pairs with free <math>\beta</math>-hCG MoAb (#70205) in the human free <math>\beta</math>-hCG sandwich immunoassay.</li><li>2. This antibody also pairs itself for in the human total hCG sandwich immunoassay.</li></ol>
<b>Note:</b>	Slight precipitation of antibody may occur on storage at 2-8°C. If necessary, refilter or centrifuge at 3,500 rpm for 15 minutes to remove the tiny precipitate.