

BIOCHECK, INC.
323 Vintage Park Drive
Foster City, CA 94404
Tel: 650-573-1968 Fax: 650-573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Polyclonal Anti-Testosterone, Protein A Purified

Catalog Number: 70656

Antigen Source: Testosterone 3-CMO-BSA

Lot Number:

Quantity: _____ mg (= _____ ml)

Concentration: _____ mg/ml

Product Form: Protein A Purified. In Borate Buffer, pH=8.0, with 0.1% NaN₃ as Preservative.

Storage: 2-8°C

Cross-Reactivity:

Testosterone	100%
Dihydro Testosterone	3.1%
Androstenedione	1.3%
EStradiol 17-b	<0.001%

Suggested Use: For Testosterone Competitive Immunoassay.