

**Product Name:** AMH Rat monoclonal antibody, clone 45A3  
**Product:** Purified Monoclonal Antibody **Catalog No.:** DM0024-B4  
**Lot Number:** **Concentration:**  
**Production Date:**  
**EXP Date:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	Elisa	Datasheet	

**Prepared By:**

**Reviewed By:**

No.168 Mashan Mellang Road | Wuxi, Jiangsu Province, China |0510-81830666