



## Certificate of Analysis

**Product Name:** TSHR Human monoclonal antibody , clone 16D10  
**Product:** Purified Monoclonal Antibody **Catalog No.:** DM0030-B2  
**Lot Number:** **Concentration:**  
**Production Date:**  
**Exp Date:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	ELISA	Datasheet	

**Prepared By:**

Quality Control-Date

**Reviewed By:**

Quality Assurance-Date

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

[www.origene.com](http://www.origene.com)

[custsupport@origene.com](mailto:custsupport@origene.com)