

**Product Name:** Sheep Anti-Thyroxin (T4) monoclonal antibody ,clone 5D2  
**Product:** Purified Monoclonal Antibody   **Catalog No.:** DM0046-B1  
**Lot Number:**                                   **Concentration:**  
**EXP Date:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	ELISA	Datasheet	

**Prepared By:** \_\_\_\_\_  
Quality Control-Date

**Reviewed By:** \_\_\_\_\_  
Quality Assurance-Date