

Product Name: Sheep Anti-Thyroxin(T4) monoclonal antibody, clone 47C7
Product: Purified Monoclonal Antibody **Catalog No.:** DM0046-B2
Lot Number: **Concentration:**
Production Date:
EXP Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Application	ELISA	Datasheet	

Prepared By:

Reviewed By:

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.com

custsupport@origene.com