

Product Name: T4 Mouse monoclonal antibody, clone 57C6
Product: Purified Monoclonal Antibody **Catalog No.:** DM0046-B3
Lot Number: **Concentration:**
Production
Expiring Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	2.77 mg/mL
Purity	SDS-Page (reduced)	G3302	>95%
Isotype	ISOSTrip™	G3346	Mouse IgG1
Application	Elisa	Datasheet	PASS

Prepared By:

Quality Control-Date

Reviewed By:

Quality/Assurance-Date

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.com custsupport@origene.com