

**Product Name:** Mouse Anti-human Abeta monoclonal antibody , clone 3G8  
**Product:** Purified Monoclonal Antibody **Catalog No.:** DM0070-B3  
**Lot Number:** **Concentration:**  
**Production Date**  
**Exp Date:**

**Antibody Purification**

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

<b>Test:</b>	<b>Method:</b>	<b>Document:</b>	<b>Results:</b>
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOStrip™	G3346	
Application	ELISA	Datasheet	

**Prepared By:****Reviewed By:**