

**Product Name:** Human Anti-human Tg monoclonal antibody , clone 12D5  
**Product:** Purified Monoclonal Antibody **Catalog No.:** DM0072-B1  
**Lot Number:** **Concentration:**  
**Production Date:**  
**Exp Date:**

**Antibody Purification**

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	ELISA	Datasheet	

**Prepared By:****Reviewed By:**