



Certificate of Analysis

Product Name: Mouse Anti-human FDP monoclonal antibody , clone 62F13
Product: Purified Monoclonal Antibody **Catalog No.:** DM0081-B1
Lot Number: **Concentration:**
Production Date
Exp Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	ELISA	Datasheet	

Prepared By:

Quality Control-Date

Reviewed By:

Quality Assurance-Date

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.com

custsupport@origene.com