



## Certificate of Analysis

**Product Name:** T4-BIGG  
**Catalog No.:** DP0042-P2  
**Lot Number:**  
**Manufactured:**  
**EXP date:**

**Concentration:**

### Antibody Purification

**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

| Test:       | Method: | Document: | Results: |
|-------------|---------|-----------|----------|
| Application | ELISA   | Datasheet |          |

**Prepared By:**

Quality Control-Date

**Reviewed By:**

Quality Assurance-Date

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

[www.origene.com](http://www.origene.com)

[custsupport@origene.com](mailto:custsupport@origene.com)