



Certificate of Analysis

Product Name: TSHR Protein
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: DP0047-P1
Concentration :
EXP date:

Antibody Purification

Buffer/Components: 0.1M Sodium citrate, pH7.2
Storage Conditions: -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page(reduced)	G3302	

Certified by:

Date:

No.168 Mashang Meiliang Road | Wuxi , Jiangsu province, China | 0510-81830666

www.origene.com

custsupport@origene.com