

## Certificate of Analysis

---

### Product Information

Product name: TPO protein

Catalog No: TP701028

Lot Number:

Quantity:

Concentration:

Buffer Formulation: PBS, pH7.4.

Storage Conditions: Store at - 80 °C. Avoid repeated freeze-thaw cycles. Stable for at least 3 years from receipt of products under proper storage and handling conditions.

---

### Quality Control

Test Items	Specifications	Results:
Protein Concentration	Microplate BCA Protein Assay	Conformed
Purity(SDS-PAGE)	> ___%	Conformed

Certified by: \_\_\_\_\_

Date: \_\_\_\_\_