



## Certificate of Analysis

**Product Name:** TG  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Manufactured:**

**Catalog No.:** TP701224-B  
**Concentration :**  
**EXP date:**

### Protein Purification

**Starting Material:** Culture supernatant  
**Buffer/Components:** PBS.  
**Storage Conditions:** -80°C

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

**Certified by:**

**Date:**

No.168 Mashang Meiliang Road | Wuxi , Jiangsu province, China | 0510-81830666

[www.origene.com](http://www.origene.com)

[custsupport@origene.com](mailto:custsupport@origene.com)